

## ADVISORY NOTICE

### **Increased risk of false positive *Proteus* results using FilmArray® Blood Culture Identification (BCID) Panel (Part No.: RFIT-ASY-0126 and RFIT-ASY-0127) with BD BACTEC™ Blood Culture Bottles - Expansion**

Dear Valued FilmArray® Customer,

The purpose of this letter is to inform you that BioFire Defense, LLC, an authorized distributor of the FilmArray® Blood Culture Identification (BCID) Panel, has been notified that BioFire Diagnostics, LLC, the legal manufacturer of the FilmArray® BCID Panel, has identified an increased risk of false positive *Proteus* results when the FilmArray BCID Panel is used with all lots of BD BACTEC™ Blood Culture Bottles (see Tables 1 & 2 below) with expiration dates of September 30, 2018 and beyond. *This is an expansion of BioFire's previous Customer Advisory Notice (FLM-PRT-0249-01) to include additional BD Blood Culture Bottle catalog numbers (see Table 2).* BioFire Diagnostics and BD are continuing to investigate the issue.


**Table 1. Previously identified media types**

<b>BD Blood Culture Bottle Catalog No.</b>	<b>Description</b>
442020	BD BACTEC™ Peds Plus™/F Culture Vials (Plastic)
442021	BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials
442023	BD BACTEC™ Plus Aerobic/F Culture Vials
442192	BD BACTEC™ Plus Aerobic/F Culture Vials
442265	BD BACTEC™ Lytic/10 Anaerobic/F Culture Vials

**Table 2. Recently identified media types**

<b>BD Blood Culture Bottle Catalog No.</b>	<b>Description</b>
442194	BD BACTEC™ Peds Plus™/F Culture Vials (Glass)
442022	BD BACTEC™ Plus Anaerobic/F Culture Vials (Plastic)
442193	BD BACTEC™ Plus Anaerobic /F Culture Vials (Glass)
442260	BD BACTEC™ Standard Aerobic/F Culture Vials (Glass)
442191	BD BACTEC™ Standard Anaerobic/F Culture Vials (Glass)

The most probable cause for this risk is the presence of an increased level of nucleic acid from non-viable *Proteus* in BD BACTEC™ Blood Culture Bottles (Tables 1 & 2). BioFire Diagnostics has confirmed the presence of *Proteus* nucleic acid in several of the lots of affected media using an independent PCR/bi-directional sequencing method. **The presence of non-viable organism does not compromise the intended function of the blood culture bottles (culturing viable microorganisms).** However, the



FilmArray BCID Panel detects nucleic acid from viable and non-viable organisms alike. Observed false positives are typically seen as multiple positives with the FilmArray BCID Panel because a positive culture is a prerequisite to a BCID test.

FilmArray BCID Panel product literature includes the following limitations: blood culture media may contain non-viable organisms and/or nucleic acid at levels that can be detected by the FilmArray BCID Panel leading to false positive test results. Typically, these false positives will be present with more than one positive result because the BCID Panel will also detect the organism that is growing in the culture bottle. In rare cases, the Gram stain result and results of the FilmArray BCID Panel may be discrepant. In these cases, the results should be used in conjunction with other clinical and laboratory findings.

**If the FilmArray BCID Panel is used to test BD BACTEC™ Blood Culture Bottles (Tables 1 & 2) with expiration dates of September 30, 2018 or later, positive results for *Proteus* should be confirmed by another method prior to reporting the test results.** BioFire Diagnostics, BioFire Defense, and BD will continue to monitor the issue in the field and will notify you if additional pertinent information is uncovered.

Please complete the accompanying *Field Acknowledgement of Receipt Form* included with this advisory notice and return via one of the listed return methods so that BioFire may acknowledge your receipt of this notification.

If you have any questions or concerns, please don't hesitate to contact our customer support department at [Support@biofiredefense.com](mailto:Support@biofiredefense.com) or via telephone by dialing 1.801.262.3592 for product technical support.

Thank you for your understanding and patience in this matter.



Robert J. Lollini  
CEO

FLM2-PRT-0253-01



79 West 4500 South, Suite 14 • Salt Lake City, Utah 84107, USA • 801.262.3592 local • 801.447.6907 fax  
[BioFireDefense.com](http://BioFireDefense.com)