



# RAZOR<sup>®</sup> Yeast Screen Assay Card

Yeast Target



**Not for Diagnostic Use**

**Kit Part No: PATH-ASY-0006**

**⚠ Handle all unknown samples and waste materials as if they were capable of transmitting infectious agents.**

This reagent kit is designed to test for DNA from **Yeast (*S. cerevisiae*)** from two unknown liquid or dry samples. These procedures are to be performed manually and have been simplified to reduce the risk of operator error.

Organism	Common Name
<i>Saccharomyces cerevisiae</i>	Yeast

Kit Contents			
5 mL Reagent Grade Water (1)	5 mL Unknown Sample Bottles (2)	1 mL Syringes with Cannula Tips (4)	Transfer Pipettes Pack (1)
Sample Swab Pack (1)	Instruction Booklet (1)	Freeze-dried Reagents in a Foil Bag (1)	

## Run Protocol Barcode

- Verify the Run Protocol is loaded onto the instrument before preparing the pouch.
- The Run Protocol may be loaded by scanning the square protocol barcode to the right into the instrument.

**NOTE:** *Times, temperatures, and results analyses are included in the protocol barcode.*



**Kit Part Number: PATH-ASY-0006**  
**Protocol Code: ITYEASTB**





# RAZOR<sup>®</sup> Yeast Screen Assay Card

Yeast Target



## Pouch Barcode

**NOTE:** *Prepare sample before loading.*

- a. After the square protocol barcode has been scanned, scan the rectangular barcode on the pouch fitment. If the rectangular barcode is damaged, use the generic one shown here.


**NOTE:** *The barcode may only be used for one run. To re-use it delete the previous run.*



ITYEASTB-ITYST001

## Sample Loading Order

- 1 Negative Port**  
Add 0.5 mL Reagent Grade Water
- 2 Unknown 1 (U1) Port**  
Add 0.5 mL Sample 1
- 3 Unknown 2 (U2) Port**  
Add 0.5 mL Sample 2
- 4 Positive Port**  
Add 0.5 mL Reagent Grade Water



## Inserting the Pouch into the RAZOR

- a. Insert the prepared pouch into the instrument with the label facing up and plungers toward the front of the instrument.

**NOTE:** *The pouch fitment should fit perfectly in the groove over the insert slot.*

