



RAZOR[®] The 10[®] (The Ten[®]) Assay Card

10 DNA Target Assay

10

Not for Diagnostic Use

Kit Part No: PATH-ASY-0061

⚠ Handle all unknown samples and waste materials as if they were capable of transmitting infectious agents.
This reagent kit is designed to test one unknown liquid or powder sample for **10 common BioThreats**.

Organism	Common Name	Organism	Common Name
<i>Bacillus anthracis</i>	Anthrax	<i>Francisella tularensis</i>	Tularemia
<i>Brucella</i> spp.	Brucella	<i>Ricinus communis</i>	Ricin
<i>Clostridium botulinum</i>	Botulism	<i>Salmonella</i> spp.	Salmonella
<i>Coxiella burnetii</i>	Q fever	<i>Variola</i>	Smallpox
<i>Escherichia coli</i> O157	E. coli O157	<i>Yersinia pestis</i>	Plague

Kit Contents			
5 mL Reagent Grade Water (1)	5 mL Unknown Sample Bottles (1)	3 mL Syringes with Cannula Tips (2)	Transfer Pipettes Pack (1)
Sample Swab Pack (1)	Instruction Booklet (1)	Freeze-dried Reagents in a Foil Bag (1)	

Run Protocol Barcode

- Verify the Run Protocol is loaded onto the instrument before preparing the pouch.
- The Run Protocol may be loaded by scanning the square protocol barcode to the right into the instrument.

NOTE: *Times, temperatures, and results analyses are included in the protocol barcode.*



Kit Part Number: PATH-ASY-0061
Protocol Code: 10TARGTB

10



RAZOR® The 10® (The Ten®) Assay Card

10 DNA Target Assay

10

Pouch Barcode

NOTE: Prepare sample before loading.

- a. After the square protocol barcode has been scanned, scan the rectangular barcode on the pouch fitment. If the rectangular barcode is damaged, use the generic one shown here.



10TARGETB-10TAR008

NOTE: The barcode may only be used for one run. To re-use it delete the previous run.


Sample Loading Order

1 Inlet Port 1

Load 2.0 mL **Unknown Sample** with 3 mL syringe

2 Inlet Port 2

Load 0.4 mL **Reagent Grade Water** from bottle with 3 mL syringe



Inserting the Pouch into the RAZOR

- a. Insert the prepared pouch into the instrument with the label facing up and plungers toward the front of the instrument.

NOTE: The pouch fitment should fit perfectly in the groove over the insert slot.

