



RAZOR[®] Water Screen Assay Card

Cryptosporidium, *E. coli* O157, and *Salmonella*



Not for Diagnostic Use

Kit Part No: PATH-ASY-0011

⚠ Handle all unknown samples and waste materials as if they were capable of transmitting infectious agents.
The Water Screen kit is designed to test liquid or dry samples for the following organisms:

Organism	Common Name	Organism	Common Name
<i>Cryptosporidium</i>	Cryptosporidium	<i>Escherichia coli</i> O157	E. coli O157
<i>Salmonella</i> spp.	Salmonella		

Kit Contents			
5 mL Reagent Grade Water (1)	5 mL Unknown Sample Bottles (2)	1 mL Syringes with Cannula Tips (4)	Transfer Pipettes Pack (1)
Sample Swab Pack (1)	Instruction Booklet (1)	Freeze-dried Reagents in a Foil Bag (1)	Loading Instructions Sticker (1)

Run Protocol Barcode

- Verify the Run Protocol is loaded onto the instrument before preparing the pouch.
- The Run Protocol may be loaded by scanning the square protocol barcode to the right into the instrument.

NOTE: *Times, temperatures, and results analyses are included in the protocol barcode.*



Kit Part Number: PATH-ASY-0011
Protocol Code: WATRSCNB





RAZOR® Water Screen Assay Card

Cryptosporidium, E. coli O157, and Salmonella



Pouch Barcode

NOTE: Prepare sample before loading.

- a. After the square protocol barcode has been scanned, scan the rectangular barcode on the pouch fitment. If the rectangular barcode is damaged, use the generic one shown here.




WATRSCNB-WATRB016

NOTE: The barcode may only be used for one run. To re-use it delete the previous run.

Sample Loading Order

- 1** Negative Port
Add 0.5 mL Reagent Grade Water
- 2** Unknown 1 (U1) Port
Add 0.5 mL Sample 1
- 3** Unknown 2 (U2) Port
Add 0.5 mL Sample 2
- 4** Positive Port
Add 0.5 mL Reagent Grade Water



Inserting the Pouch into the RAZOR

- a. Insert the prepared pouch into the instrument with the label facing up and plungers toward the front of the instrument.

NOTE: The pouch fitment should fit perfectly in the groove over the insert slot.

