



Book No. ASAY-PRT-0131

# *Freeze-Dried Reagent*

## *Detection Kit Instruction Booklet*

For Hybridization Probe Assays



Manufactured by  
BioFire Defense, LLC  
79 West 4500 South, Suite 14  
Salt Lake City, Utah 84107 USA  
+1-801-262-3592 | [www.BioFireDefense.com](http://www.BioFireDefense.com)

## CUSTOMER AND TECHNICAL SUPPORT

### Reach Us on the Web

<http://www.BioFireDefense.com>

### Reach Us by E-mail

[support@biofiredefense.com](mailto:support@biofiredefense.com)

### Reach Us by Phone

8 a.m. to 5 p.m. - Mt. Standard Time

1-801-262-3592 - US and Canada  
IDD Prefix-1-801-262-3592 - Intl

### Reach Us by Fax

1-801-967-3688 - US and Canada  
IDD Prefix-1-801-967-3688 - Intl

© Copyright 2006-2017, BioFire Defense  
All rights reserved.

Hybridization Probe Freeze-Dried Reagent Instruction Booklet

ASAY-PRT-0131-03 01/17

The information contained in this document is subject to change without notice. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of BioFire Defense, LLC.

Platinum Path and IT 1-2-3 are trademarks of BioFire Defense, LLC or BioFire Diagnostics, LLC. The BioFire logo is a registered trademark of BioFire Defense LLC or BioFire Diagnostics, LLC. All other names of products and brands appearing in this manual are trademarks or registered trademarks of their respective owners.

The purchase of this product includes a limited, nontransferable license under specific claims of one or more U.S. patents as listed on BioFire Defense's web site (<http://www.BioFireDefense.com/LegalNotices/>) and owned or licensed by the University of Utah Research Foundation and/or BioFire, to use only the enclosed amount of product according to the specified protocols. No right is conveyed, expressly, by implication, or by estoppel, to use any instrument or system under any claim of such U.S. patent(s), other than for the amount of product contained herein.

The purchase of this product allows the purchaser to use it for amplification and detection of nucleic acid sequences for providing human in vitro diagnostics. No general patent or license of any kind other than this specific right of use from purchase is granted hereby.

BioFire Defense, LLC makes no warranty of any kind with regard to this material, including, but not limited to the implied warranties of merchantability and fitness for a particular purpose. BioFire Defense, LLC shall not be liable for errors contained herein or for incidental consequential damages in connection with the furnishing, performance or use of this material.

Components of ASAY-ASY-0210 are licensed from MO BIO Laboratories, Inc. under certain patents owned or controlled by MO BIO Laboratories, Inc. for use in ASAY-ASY-0120. No other rights, express or implied, are conveyed by the sale of this product.

This product is for presumptive testing only, and is not to be used for diagnosis of disease. Positive tests should be confirmed by a certified public health laboratory. Negative results are possible even when the target organism is present, due to various factors including PCR inhibitors and improper sample preparation, damage to the test kit or instrument, or human error.

# Table of Contents

Introduction: Kit Contents.....	2
General Safety Precautions .....	2
Laboratory Precautions.....	3
Storage .....	3
Preventing Contamination During Setup .....	3
Product Overview.....	4
Basic Stages of the Protocol .....	4
Reagent Setup.....	5
Negative Controls .....	5
Unknown Samples.....	6
Positive Controls.....	7
Centrifuging Capillaries.....	7
Amplification Procedures .....	8
Two Types of Tests .....	8
R.A.P.I.D. 7200 Protocol (DNA).....	9
R.A.P.I.D. 7200 Protocol (RNA).....	11
R.A.P.I.D. 9200 and LT Protocol (DNA).....	13
R.A.P.I.D. 9200 and LT Protocol (RNA).....	14
Results .....	16
Appendix A: Decontamination.....	19
Broken Capillary .....	19
Sample Carousel Decontamination .....	20
Instrument Decontamination.....	21
Decontamination for Other Areas .....	21
Appendix B: Troubleshooting .....	22
Appendix C: Melt Protocol .....	23
R.A.P.I.D. 7200 Melt Protocol.....	23
R.A.P.I.D. 9200 and LT Melt Protocol.....	23
Appendix D: Ordering Information .....	24
Freeze-Dried Reagent Kits .....	24
Sample Purification Kits and Supplies.....	25



## Introduction: Kit Contents

These reagents are for the identification of their described target agent. Use the accompanying assay card for specific information relating to the described target.

This instruction booklet describes procedures to perform real-time polymerase chain reaction (PCR) based assays on BioFire Defense's Ruggedized Advanced Pathogen Identification Device (R.A.P.I.D.®) Instruments. These procedures describe reagent setup and results analysis. For additional information, consult the system manual that accompanied the instrument and the accompanying assay card.

<b>Kit Contents</b>		
<b>Vial</b>	<b>Label</b>	<b>Contents and Use</b>
<b>Green Label</b>	Positive Control	Ready-to-use hybridization probe mix. Each vial contains reagents for two 20 µL reactions.  Contains purified amplicon for positive result.
<b>Red Label</b>	Negative Control	Ready-to-use hybridization probe mix. Each vial contains reagents for two 20 µL reactions.
<b>Blue Label</b>	Unknown	Ready-to-use hybridization probe mix. Each vial contains reagents for two 20 µL reactions.
<b>Amber Tube</b>	Reagent grade water	2 x 850 µL

<b>Additional Equipment Needed</b>	
R.A.P.I.D. Instrument with laptop and current software	LightCycler® capillaries

## General Safety Precautions

One of the most important rules of working in a laboratory environment is to avoid contamination: 1) use filter tips on your pipettes when working with any liquid solution, and 2) use different pipettes for setting up the reactions and for handling the samples.



## Laboratory Precautions

For general biosafety guidelines, refer to *Biosafety in Microbiological and Biomedical Laboratories* (BMBL), 4th Edition, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention and National Institutes of Health (DHHS), May 1999. Available from <http://www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4toc.htm>

### Handling of Biohazard Wastes

Dispose of used reagent vials in accordance with good laboratory practices. Waste from possible positive samples, or samples containing amplified DNA molecules, should be inactivated using appropriate procedures before disposal.

For more information refer to DHHS National Institute of Occupational Safety and Health (NIOSH) Publication No. 88-119, *Guidelines for Protecting the Safety and Health of Health Care Workers* (section 6, “Hazardous Waste Disposal”).

### Storage

Store contents in a dry environment at room temperature (18–30 °C). **DO NOT REFRIGERATE.** Do not allow the kit to sit in direct sunlight. When pouches are opened, use contents immediately. Discard unused reagents from open pouches. Reagents are stable through the expiration date printed on the label provided the storage recommendations have been followed. **Do not use reagents after they have expired.**

### Preventing Contamination During Setup

Some concerns for users of PCR based assays for detection are the sensitivity levels of the assay and the high copy number of product that is generated from the amplification. Both of these factors raise the risk of false positive calls, especially if a capillary is broken. To mitigate this risk of contamination, some important precautions can be taken. These precautions include collecting, purifying, preparing, and running a test in separate rooms or areas, as well as executing decontamination procedures after a capillary breakage occurs.



Outside Sample Collection  
of Raw Sample

Lab Room / Area 1  
Sample Purification

Lab Room / Area 2  
Reaction Prep

Lab Room / Area 3  
Run Reaction



Start

Finish

## Product Overview

This kit is specially adapted for amplification in glass capillaries using the R.A.P.I.D. Instrument and hybridization probes.

A fragment of the target DNA is amplified using specific primers. The amplicon is detected by fluorescence using a specific pair of hybridization probes. These probes consist of two different short oligonucleotides that hybridize to an internal sequence of the amplified fragment during the annealing phase of the reaction cycle. One probe is labeled at the 5' end with LCRed 640. To avoid extension on the 3' end, it is modified by phosphorylation. The second probe is labeled at the 3' end with fluorescein. Only after hybridization to the template DNA do the two probes come in close proximity, resulting in fluorescence resonance energy transfer (FRET) between the two fluorophores. The fluorescence emitted by the LCRed 640 dye is measured in channel 2 of the R.A.P.I.D. Instrument.

The fluorescent signal from the unknown sample is compared to the signals from the positive and negative control samples.


### Basic Stages of the Protocol


Amplify and monitor online in real-time for the presence of the target DNA.

**Step 1:** Purify and mix DNA with freeze-dried reagents and water

**Step 2:** Test samples using the R.A.P.I.D. Instrument

**Step 3:** View results with the auto-call software



 **Note:** BioFire Defense strongly recommends using established DNA purification procedures using their IT 1-2-3™ DNA Purification Kits and a Vortex-Genie® 2 with 2 mL Turbomix® adapter.

### Assay Time

Based on an assay for 2 positive, 2 negative, and 28 unknown samples:

Procedure	Time
Set up test	15 min.
Run test using R.A.P.I.D. Instrument	25 min.
<b>Total Experiment Time</b>	<b>40 min.</b>

### Capillary Handling

Never touch the reaction capillaries with bare hands. Always wear gloves when handling capillaries.

### Reagent Setup

Before preparing reagents, obtain the appropriate numbers of Negative Control (NC), Positive Control (PC), and Unknown reagent vials, as well as Reagent Grade Water. Each reagent vial contains enough master mix for two reactions. Freeze-dried reagent vials contain all of the components required for PCR. Negative and Unknown vials lack DNA template and are used to setup NC and Unknown samples. PC vials contain reagents that have been freeze-dried with specified quantities of target template DNA.

Reactions should be prepared so that the NC samples are set up first, then the Unknown samples, and lastly the PC samples. Preparation in this order will reduce the potential of contaminating the Unknown samples with PC template.

### Negative Controls

Put on a clean pair of gloves. Prepare the freeze-dried NC samples as follows:

1. To ensure that the reagent pellet or pellet fragments are at or near the bottom of the vial, tap the NC reagent vial on the bench or centrifuge for 3 sec.
2. Remove the rubber cap from each NC reagent vial.




3. Visually check that the reagent pellet is at or near the bottom of the vial. If not, with a clean pipette tip, push the pellet or fragment(s) to the bottom of the vial. Change the pipette tip.
4. To the NC reagent vial, add **40  $\mu\text{L}$  of Reagent Grade Water**. Change the pipette tip.
5. Replace cap on each reagent vial.
6. Vortex the NC reagent vial for 5 sec at maximum speed.
7. Load reagent vial in minicentrifuge and centrifuge for 3 sec to bring all liquid to the bottom of the vial.
8. Visually check that the entire pellet has been rehydrated. If not repeat steps 6 through 8.
9. Pipette **19  $\mu\text{L}$  of the hydrated mixture** into a capillary. Repeat pipette step into a second capillary to obtain a duplicate.
10. Cap capillaries and discard gloves. Record capillary number, sample type, and assay name. Change the pipette tip.

## Unknown Samples

---

Put on a clean pair of gloves. Prepare freeze-dried Unknown samples as follows:

 **Note:** To minimize the risk of cross-contamination, each unknown sample should be processed from steps 1 through 11 before another unknown sample is set up.

1. To ensure that the reagent pellet or pellet fragments are at or near the bottom of the vial, tap the Unknown reagent vial on the bench or centrifuge for 3 sec.
2. Remove the rubber cap from the Unknown reagent vial.
3. Visually check that the reagent pellet is at or near the bottom of the vial. If not, with a clean pipette tip, push the pellet or fragment(s) to the bottom of the vial. Change the pipette tip.
4. To the Unknown reagent vial, add **40  $\mu\text{L}$  of Sample**.
5. Replace cap on reagent vial.
6. Vortex the Unknown reagent vial for 5 sec at maximum speed.
7. Load Unknown reagent vial in minicentrifuge and centrifuge for 3 sec to bring all liquid to the bottom of the vial.
8. Visually check that the entire pellet has been rehydrated. If not repeat steps 6 through 8.
9. Pipette **19  $\mu\text{L}$  of the hydrated mixture** into a capillary. Repeat pipette step into a second capillary to obtain a duplicate. Change the pipette tip.





11. Cap capillaries and discard gloves. Record capillary number, sample information, and assay name.

## Positive Controls

---

**⚠ WARNING:** High risk of cross-contamination. Follow instructions carefully.

Put on a clean pair of gloves. Prepare freeze-dried PC samples as follows:

1. To ensure that the reagent pellet or pellet fragments are at or near the bottom of the vial, tap the PC reagent vial on the bench or centrifuge for 3 sec.
2. Remove the rubber cap from each PC reagent vial.
3. Visually check that the reagent pellet is at or near the bottom of the vial. If not, with a clean pipette tip, push the pellet or fragment(s) to the bottom of the vial. Change the pipette tip.
4. To the PC reagent vial, add **40  $\mu$ L of Reagent Grade Water**. Change the pipette tip.
5. Replace cap on reagent vial.
6. Vortex the PC reagent vial for 5 sec at maximum speed.
7. Load reagent vial in minicentrifuge and centrifuge for 3 sec to bring all liquid to the bottom of the vial.
10. Visually check that the entire pellet has been rehydrated. If not repeat steps 6 through 8.
11. Pipette **19  $\mu$ L** of the **hydrated mixture** into a capillary. Repeat this step into a second capillary to obtain a duplicate.
12. Cap capillaries and discard gloves. Record capillary number, sample type, and assay name. Discard pipette tip.

## Centrifuging Capillaries

---

1. Once all of the reagents have been prepared and transferred to capillaries, load capillaries in a minicentrifuge fitted with capillary adaptors.
2. Centrifuge at the **lowest speed setting** by holding down the pulse button for 3 sec. This transfers the liquid down to the tips of the capillaries and prepares them for testing.

**⚠ WARNING:** Centrifuging capillaries at full speed will cause them to break.



3. Once the samples are prepared, use the following protocol for your specific instrument. Please note especially that the denature temperatures are different between the R.A.P.I.D. 7200 and R.A.P.I.D. 9200. Follow the instructions found in the instrument users guide on loading the carousel and running a test.

## Amplification Procedures

---

### Starting the R.A.P.I.D. Software

Refer to the R.A.P.I.D. System Manual for detailed instructions on setting up and using the R.A.P.I.D. software for an experiment.

### Experimental Protocol

This protocol consists of three programs. Use the appropriate protocol for your instrument.

**Program 1:** Denaturation of the target DNA

**Program 2:** Amplification of the target DNA

**Program 3:** Melting of the target DNA

## Two Types of Tests

---

**Batch Tests**—Batch tests are used for a single assay on multiple samples. Follow the directions below to avoid cross-contamination.

**Screen Tests**—Screen tests are for multiple assays on a single or limited number of samples. To minimize the risk of cross-contamination, prepare all Unknown reactions from one purified sample and transfer to capillaries before beginning setup with another unknown sample. Next, prepare all NCs and transfer them to capillaries. Prepare all PCs last and transfer to capillaries to finish the reconstitution procedure.



## R.A.P.I.D. 7200 Protocol (DNA)

The following protocol is for the R.A.P.I.D. 7200 Instrument. The image below illustrates where the protocol needs to be entered in the run software.

The screenshot shows the software interface for entering protocol data. It is divided into two main sections: 'Cycle Program Data' and 'Temperature Targets'.  
**Cycle Program Data:** 'Cycles' is set to 1, and 'Analysis Mode' is set to None.  
**Temperature Targets:** This section lists several parameters with corresponding input fields: Target Temperature (94 °C), Incubation Time (2:00 hrs:min:sec), Temperature Transition Rate (20.0 °C/sec), Secondary Target Temperature (0 °C), Step Size (0 °C), Step Delay (0 cycles), and Acquisition Mode (NONE). There are also 'Add' and 'Del' buttons at the bottom of the input fields.

### Program 1: DNA Denaturation

Cycle Program Data		Value
Cycles		1
Analysis Mode		None
Temperature Targets		Segment 1
Target Temperature		94 °C
Incubation Time (hrs:min:s)		120 s
Temperature Transition Rate		20 °C/s
Second Target Temperature		0 °C
Step Size		0 °C
Step Delay (Cycles)		0
Acquisition Mode		None



**Temperature Targets**

Target Temperature (°C)  
 Incubation Time (hrs:min:sec)  
 Temperature Transition Rate (°C / sec)  
 Secondary Target Temperature (°C)  
 Step Size (°C)  
 Step Delay (cycles)  
 Acquisition Mode

94 0 20.0 0 0.0 0 NONE   
 60 20 20.0 0 0.0 0 SINGLE

## Program 2: Amplification


Cycle Program Data		Value	
Cycles	45		
Analysis Mode	Auto Analysis		
Temperature Targets	Segment 1	Segment 2	
Target Temperature	94 °C	60 °C	
Incubation Time (hrs:min:s)	0 s	20 s	
Temperature Transition Rate	20 °C/s	20 °C/s	
Second Target Temperature	0 °C	0 °C	
Step Size	0 °C	0 °C	
Step Delay (Cycles)	0	0	
Acquisition Mode	None	Single	

## Fluorescence Parameters

**Display Mode:** fluorescence channel 2 (Ch2/1)

Set the fluorescence gains as follows:

Fluorimeter Gain	Value
Channel 1 (F1)	1
Channel 2 (F2)	8
Channel 3 (F3)	1



## R.A.P.I.D. 7200 Protocol (RNA)

**Cycle Program Data**

Cycles:

Analysis Mode:

---

**Temperature Targets**

Target Temperature (°C)

Incubation Time (hrs:min:sec)

Temperature Transition Rate (°C/sec)

Secondary Target Temperature (°C)

Step Size (°C)

Step Delay (cycles)

Acquisition Mode

40  20.0  0.0  NONE

### Program 1: RNA Sample Preparation Hold

Cycle Program Data		Value
Cycles		1
Analysis Mode		None
Temperature Targets		Segment 1
Target Temperature		40 °C
Incubation Time (hrs:min:s)		30 min
Temperature Transition Rate		20 °C/s
Second Target Temperature		0 °C
Step Size		0 °C
Step Delay (Cycles)		0
Acquisition Mode		None



### Program 2: RNA Denaturation

Cycle Program Data		Value
Cycles		1
Analysis Mode		None
Temperature Targets		Segment 1
Target Temperature		94 °C
Incubation Time (hrs:min:s)		120 s
Temperature Transition Rate		20 °C/s
Second Target Temperature		0 °C
Step Size		0 °C
Step Delay (Cycles)		0
Acquisition Mode		None

### Program 3: RNA Amplification

Cycle Program Data		Value	
Cycles		45	
Analysis Mode		Auto Analysis	
Temperature Targets		Segment 1	Segment 2
Target Temperature		94 °C	60 °C
Incubation Time (hrs:min:s)		0 s	20 s
Temperature Transition Rate		20 °C/s	20 °C/s
Second Target Temperature		0 °C	0 °C
Step Size		0 °C	0 °C
Step Delay (Cycles)		0	0
Acquisition Mode		None	Single



## Fluorescence Parameters

**Display Mode:** fluorescence channel 2 (Ch2/1)

Set the fluorescence gains as follows:

Fluorimeter Gain	Value
Channel 1 (F1)	1
Channel 2 (F2)	8
Channel 3 (F3)	1

## R.A.P.I.D. 9200 and LT Protocol (DNA)

The following protocol steps are for the R.A.P.I.D. 9200 and LT Instruments. The image below illustrates where the protocol needs to be entered.

The screenshot shows the software interface for setting a protocol. At the top, there is a dropdown menu for 'Color Comp (Off)' set to 'RP\_3 Ch. on CDM1 (RP\_)'. Below this is a 'Notes' section with a 'Programs' table. The 'Programs' table has columns for 'Program Name', 'Cycles', and 'Analysis Mode'. One program is listed with 'Program' as the name, '1' for cycles, and 'None' for analysis mode. Below the programs section is a 'Temperature Targets' section with a table. The 'Temperature Targets' table has columns for 'Target (°C)', 'Hold (hh:mm:ss)', 'Ramp Rate (°C/s)', 'Sec. Target (°C)', 'Step Size (°C)', 'Step Delay (cycles)', and 'Acquisition Mode'. One target is listed with '92' for target, '00:02:00' for hold, '20' for ramp rate, and '0' for the other three columns, with 'None' for acquisition mode.

### Program 1: DNA Denaturation

Cycle Program Data	Value
Cycles	1
Analysis Mode	None
Temperature Targets	Segment 1
Target	92 °C
Hold (hh:mm:ss)	120 s
Ramp Rate	20 °C/s



Sec Target	0 °C
Step Size	0 °C
Step Delay (cycles)	0
Acquisition Mode	None

## Program 2: DNA Amplification

Cycle Program Data		Value	
Cycles	45		
Analysis Mode	Quantification		
Temperature Targets		Segment 1	Segment 2
Target	92 °C	60 °C	
Hold (hh:mm:ss)	0 s	20 s	
Ramp Rate	20 °C/s	20 °C/s	
Sec Target	0 °C	0 °C	
Step Size	0 °C	0 °C	
Step Delay (cycles)	0	0	
Acquisition Mode	None	Single	


## R.A.P.I.D. 9200 and LT Protocol (RNA)

The screenshot shows the software interface for the R.A.P.I.D. 9200. At the top, the 'Editor Comp (0/0)' dropdown is set to 'RP\_3 Ch. on CDM1 (RP\_)'. Below this, the 'in Notes' section is visible. The main area contains two tables:

Programs							
Program Name	Cycles	Analysis Mode					
Program	1	None					

Temperature Targets							
Target (°C)	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Sec Target (°C)	Step Size (°C)	Step Delay (cycles)	Acquisition Mode	
40	00:30:00	20	0	0	0	None	






### Program 1: RNA Sample Preparation Hold

Cycle Program Data		Value
Cycles		1
Analysis Mode		None
Temperature Targets		Segment 1
Target Temperature		40 °C
Incubation Time (hrs:min:s)		30 min
Temperature Transition Rate		20 °C/s
Second Target Temperature		0 °C
Step Size		0 °C
Step Delay (Cycles)		0
Acquisition Mode		None

### Program 2: RNA Denaturation

Cycle Program Data		Value
Cycles		1
Analysis Mode		None
Temperature Targets		Segment 1
Target Temperature		92 °C
Incubation Time (hrs:min:s)		120 s
Temperature Transition Rate		20 °C/s
Second Target Temperature		0 °C
Step Size		0 °C
Step Delay (Cycles)		0
Acquisition Mode		None



### Program 3: RNA Amplification

Cycle Program Data		Value	
Cycles	45		
Analysis Mode	Auto Analysis		
Temperature Targets		Segment 1	Segment 2
Target Temperature	92 °C	60 °C	
Incubation Time (hrs:min:s)	0 s	20 s	
Temperature Transition Rate	20 °C/s	20 °C/s	
Second Target Temperature	0 °C	0 °C	
Step Size	0 °C	0 °C	
Step Delay (Cycles)	0	0	
Acquisition Mode	None	Single	

## Results

The results shown below are for R.A.P.I.D. Instrument models 7200, 9200, and LT. Use the appropriate instructions for your system.

### R.A.P.I.D. 7200

A successful test yields results similar to the ones shown below. Results are displayed in grouped panels if multiple organisms were tested.

Organism Present				
#	Sample Name	Control	Score	Result
1	Sample 1	Positive	1148	Present
2	Sample 2	Positive	1046	Present
4	Sample 4	Negative	1	Not Detected
5	Sample 5	Negative	2	Not Detected
6	Sample 6	Unknown	214	Present
7	Sample 7	Unknown	15	Not Detected
8	Sample 8	Unknown	6	Not Detected
9	Sample 9	Unknown	4	Not Detected



For each capillary, the results include the following:

- **#**—Identifies the carousel position of the capillary.
- **Sample Name**—Provides the name of the sample.
- **Control**—Identifies the sample type for each sample.
- **Score**—The score is derived by calculating the relationship between the sample's fluorescence and the level of background fluorescence in the sample.
- **Result**—Displays the overall result for the sample.

The possible results include the following:

<b>Present</b>	A red <b>Present</b> call for an unknown indicates that the target was identified in that sample and the controls were successful. A blue Present call is for a positive control.
<b>Not Detected</b>	A green <b>Not Detected</b> call for an Unknown indicates that the target was not identified in that sample and the test controls were successful. A blue Not Detected call is for a negative control.
<b>Please Repeat</b>	A <b>Please Repeat</b> indicates that the positive or negative control failed.


**Positive** and **Please Repeat** results should be referred to a supervisor or advanced operator for appropriate follow-up.

To display the real-time fluorescence data for each sample, select the sample of interest and click the **Show Graph** button at the bottom of the screen. For more detailed information on performing a test and analyzing the results, see the *R.A.P.I.D. System Manual* that came with your instrument.

#### **R.A.P.I.D. 9200 and LT**

A successful test yields results similar to the ones shown below. Results are displayed in grouped panels if multiple organisms were tested.



R.A.P.I.D.		
 <h2 style="display: inline;">Analysis Results</h2>		
Accession No	Test	Combined
123	B. anthracis Target 2	Negative
123	Y. pestis Target 1	Invalid
123	B. anthracis Target 3	Positive
123	B. anthracis Target 1	Invalid

For each sample, the results include the following:

- **Accession No.**
- **Test**—pathogen target
- **Combined Result**—is a combination of the control results and the unknown results. See chart below
- **Concentrations**—quantification results for certain positive samples (specific target/specimen combinations) including appropriate units.

The Results window displays results only for the capillaries that contain unknown samples.


The possible results include the following:

<b>Positive</b>	A red positive call indicates that the target was identified in that sample and the run controls were successful. A sample that is detected as positive will generate an audible signal in the software.
<b>Negative</b>	A green negative call indicates that the target was not identified in that sample and all the run controls (NC and PC) were successful.
<b>Uncertain</b>	A yellow uncertain call indicates that the instrument was unable to interpret data and the controls were successful.
<b>Invalid</b>	A yellow invalid call indicates that a PC or NC failed. All samples connected to that control would be invalid.

**Positive**, **Invalid**, and **Uncertain** results should be referred to a supervisor or advanced operator for appropriate follow-up.



Click **Finish** to view a report of the test results. For more detailed information on performing a test and analyzing the results, see the *R.A.P.I.D. System Manual* that came with your instrument.

 **Note:** Melting curve temperatures can be affected by final DNA concentration as well as contaminants like residual ethanol from DNA preparation. Temperatures are approximate and can shift depending on experimental variables.


## Appendix A: Decontamination

The list below provides items that are necessary in a laboratory for decontamination purposes.

- 10% bleach solution in a squeeze or spray bottle
- Water in a squeeze or spray bottle
- 95% ethanol in a squeeze bottle
- DNAZap™ or equivalent DNA degrading system
- Paper towels
- Nylon brush
- Bleach wipes
- Lens cleaner
- Lens paper

### Broken Capillary

---

 **WARNING:** Immediately after capillary breakage make sure no one uses any potentially contaminated areas or instruments. Follow the steps below to properly decontaminate and dispose of the capillary.

1. Put on PPE such as gloves and safety glasses.
2. Stop the run and unplug the instrument before beginning cleanup.
3. Ensure no one else uses any potentially contaminated areas or instruments.
4. Decontaminate and dispose of the broken capillary using the following steps:
  - a. Dispose of potentially contaminated gloves.
  - b. Put on a clean pair of gloves and dispose of the potentially contaminated lab coat.



- c. Clean up broken glass, discard in capillary containment bin or biohazardous waste container, and complete removal of capillaries from carousel (if necessary).
- d. Change gloves.
- e. Follow decontamination procedures for cleaning affected area.

## Sample Carousel Decontamination


---

R.A.P.I.D. 7200 carousels are custom fitted to each instrument. R.A.P.I.D. 9200 and LT carousels are interchangeable and disposable. Contact BioFire Defense to obtain a replacement carousel if the decontamination procedures fail to eliminate all contamination or the carousel is damaged.

1. Place carousel on bench top.
2. Spray 10% bleach solution on top and bottom of carousel. Let sit for 5–10 min.
3. Wipe carousel and bench area with clean paper towel.
4. Repeat steps 2 and 3 twice, for a total of three wipes.
5. Wipe carousel and bench area with new paper towel.
6. Spray distilled water on both sides of carousel.
7. Wipe carousel and bench area; dry with a new paper towel.
8. Spray both sides of carousel with DNAzap, or equivalent product. See product for proper usage.
9. Rinse by spraying distilled water on both sides of carousel.
10. Wipe carousel and bench area dry with a new paper towel. Tap the carousel on your hand or clean bench top to remove excess distilled water.
11. Use a squirt bottle of 95% ethanol to cover the entire carousel with ethanol. In particular, squirt ethanol into the capillary holes to remove any excess bleach.
12. Wipe carousel and bench area dry with a new paper towel. Tap the carousel on your hand or clean bench top to remove excess ethanol.
13. Use a new nylon brush to dry out all holes in the carousel. Make sure carousel is dry, as residual bleach could alter fluorescence readings.\* If desired, check all carousel holes by inserting a clean, empty capillary. If insertion and removal of a capillary is difficult for any hole, repeat rinsing and drying steps.
14. Discard nylon brush.

---

\* Due to brass components, any bleach that is not completely rinsed away can cause corrosion. This in turn can cause increased friction leading to broken capillaries and contamination in future runs. BioFire Defense recommends replacing disposable carousels (models 9200 and LT). Custom-fitted carousels (model 7200) must be thoroughly rinsed and dried.



## Instrument Decontamination

---

1. Remove carousel from machine and follow carousel decontamination procedures.
2. Remove broken glass from instrument (including the carousel chamber) with lens paper lightly moistened with water. You may carefully remove glass from the optics with the lens paper.
3. To remove glass from under the chamber fan, remove the fan blade according to the instrument manual.
4. Change gloves and wet a paper towel with 10% bleach and wipe the entire instrument (including the bottom of the instrument). Do not wipe bleach on the optics.
5. Repeat bleach step with a fresh paper towel and clean gloves.
6. Wet a paper towel with water and wipe the entire instrument (including the bottom of the instrument). Do not wipe the optics.
7. Repeat water step with a fresh paper towel.
8. Clean the optics of the instrument with lens cleaner and lens paper. Use light pressure when cleaning the optics.
9. Repeat optics cleaning step with a fresh lens paper.
10. Make a visual check of the thermocouple to make sure it is straight and not bent to one side or the other. This will ensure the thermocouple will not break the next carousel of capillaries. In the 7200, you can check the thermocouple alignment by placing the white storage carousel in the carousel chamber, making sure to align the notch in the storage carousel with the thermocouple.
11. Dispose of gloves and other potentially contaminated items.

## Decontamination for Other Areas

---

1. Spray 10% bleach solution on area of potential contamination. Let stand for 5 min.
2. Wipe with paper towel.
3. Repeat steps 1 and 2 twice, for a total of three wipes.
4. Spray area with distilled water.
5. Wipe dry with new paper towel.
6. Spray area with DNAZap, or equivalent product. See product for proper usage.
7. Rinse by spraying area with distilled water and wiping dry.



## Appendix B: Troubleshooting

**Problem:** *No amplifications can be monitored.*

- Possible Cause**—Wrong channel chosen for real-time monitoring  
**Recommendation**—Check that the fluorescence channel selected is Channel 2.
- Possible Cause**—Measurements do not occur.  
**Recommendation**—Make sure “single” is selected as acquisition mode in programming screen.
- Possible Cause**—Inhibitory effects of sample material.  
**Recommendation**—Purify DNA sample using recommended kits.

**Problem:** *Fluorescence intensity too high.*


**Possible Cause**—Unsuitable gain settings.  
**Recommendation**—Check gain settings in the programming screen. Lower settings if necessary.

**Problem:** *Fluorescence intensity too low.*

- Possible Cause**—Unsuitable gain settings.  
**Recommendation**—Check gain settings in the programming screen. Raise settings if necessary.
- Possible Cause**—Reaction mix exposed to light.  
**Recommendation**—Be sure to protect reaction vials from light. Try a new batch of reagents.
- Possible Cause**—Very low starting DNA.  
**Recommendation**—Use more sample DNA.

**Problem:** *Negative control samples are positive.*

**Possible Cause**—Contamination.  
**Recommendation**—Replace all critical solutions.  
Decontaminate pipettors.  
Use aerosol free tips.  
Use a fresh tip for each pipetting step.





## Appendix C: Melt Protocol

The following melt protocols are used to confirm test results.

### R.A.P.I.D. 7200 Melt Protocol

Cycle Program Data		Value		
Cycles	1			
Analysis Mode	Melting Curves			
Temperature Targets	Segment 1	Segment 2	Segment 3	
Target Temperature	94 °C	45 °C	80 °C	
Incubation Time (hrs:min:s)	2 s	140 s	0 s	
Temperature Transition Rate	20 °C/s	20 °C/s	0.2 °C/s	
Second Target Temperature	0 °C	0 °C	0 °C	
Step Size	0 °C	0 °C	0 °C	
Step Delay (Cycles)	0	0	0	
Acquisition Mode	None	None	Step	

### R.A.P.I.D. 9200 and LT Melt Protocol

Cycle Program Data		Value		
Cycles	1			
Analysis Mode	Melting Curves			
Temperature Targets	Segment 1	Segment 2	Segment 3	
Target Temperature	92 °C	45 °C	80 °C	
Incubation Time (hrs:min:s)	2 s	140 s	0 s	
Temperature Transition Rate	20 °C/s	20 °C/s	0.2 °C/s	
Second Target Temperature	0 °C	0 °C	0 °C	
Step Size	0 °C	0 °C	0 °C	
Step Delay (Cycles)	0	0	0	
Acquisition Mode	None	None	Step	



## Fluorescence Parameters

**Display Mode:** fluorescence channel 2 (Ch2/1)

Set the fluorescence gains as follows:

Fluorimeter Gain	Value
Channel 1 (F1)	1
Channel 2 (F2)	8
Channel 3 (F3)	1

## Appendix D: Ordering Information

For additional product information visit us online at [www.BioFireDefense.com](http://www.BioFireDefense.com).

### Freeze-Dried Reagent Kits

Item	Part Number
BioThreat Screening Kit Tests for: <i>B. anthracis</i> (Target 1), <i>F. tularensis</i> (Target 1), <i>Y. pestis</i> (Target 1), <i>Brucella species</i> (Target 1)	3833
Pathogen test Kit Tests for: <i>Listeria monocytogenes</i> (Target 1), <i>E. coli</i> O157 (Target 1), <i>Salmonella species</i> (Target 1), <i>Campylobacter species</i> (Target 1)	3834
Reagent Training Kit Freeze-dried reagents used for training on the R.A.P.I.D. instrument.	3832
<i>B. anthracis</i> (Anthrax, Target 1)	3828
<i>B. anthracis</i> (Anthrax, Target 2)	3841
<i>B. anthracis</i> (Anthrax, Target 3)	3887
<i>F. tularensis</i> (Tularemia, Target 1)	3830
<i>F. tularensis</i> (Tularemia, Target 2)	3882
<i>Y. pestis</i> (Plague, Target 1)	3831
<i>Y. pestis</i> (Plague, Target 2)	3881
<i>Brucella species</i> (Brucellosis, Target 1)	3829
<i>Listeria monocytogenes</i> (Listeriosis, Target 1)	3827
<i>E. coli</i> O157 (Target 1)	3825
<i>Salmonella species</i> (Salmonellosis, Target 1)	3826



Item	Part Number
<i>Campylobacter</i> species (Target 1)	3839
<i>C. botulinum</i> Type A (Botulism, Target 1)	3840
<i>Cryptosporidium</i> (Target 2)	ASAY-ASY-0101
Variola (Small Pox)	ASAY-ASY-0102
Ricin Target 1	ASAY-ASY-0103
Ricin Target 2	ASAY-ASY-0104
Avian Influenza H5 Subtype Target 1	ASAY-ASY-0105
Avian Influenza H5 Subtype Target 2	ASAY-ASY-0106
Influenza A Target 1	ASAY-ASY-0109

## Sample Purification Kits and Supplies

Item	Contents	Part Number
IT 1-2-3™ PLATINUM Path Sample Purification Kit	Supplies to purify 40 samples from a wide array of samples.	ASAY-ASY-0120
IT 1-2-3™ SWIPE Sample Purification Kit	Supplies to purify 40 samples, 18 reactions per sample	ASAY-ASY-0005
IT 1-2-3™ FLOW Sample Purification Kit	Supplies to purify 40 samples, 70 reactions per sample, includes protease	ASAY-ASY-0004
IT 1-2-3™ SCOOP Sample Purification Kit	Supplies to purify 50 samples, 18 reactions per sample	ASAY-ASY-0502
ITI 1-2-3™ VIBE Sample Purification Kit	Supplies to purify 40 samples, 9–18 reactions per sample, includes protease	ASAY-ASY-0500
ITI 1-2-3™ RNA Module	Carrier RNA	ASAY-ASY-0501
Filtered Blender Bags	1 Bag of 10 (Brinkmann)	ASAY-ASY-0060
PBS packets	Phosphate Buffered Saline, pH 7.4 PBS powder, 10 packets per box (each packet makes 1 L PBS) (Sigma-Aldrich)	ASAY-ASY-0061
Triton X-100	100 mL Triton (Sigma-Aldrich)	ASAY-ASY-0062
Stool Collection Vials (Cary Blair)	Box of 20 (PROTOCOL-Fisher)	ASAY-ASY-0063
Stool Concentrator Filters	Box of 30 (PARAPAK, Meridian Diagnostics)	ASAY-ASY-0064





79 West 4500 South, Suite 14, Salt Lake City, Utah 84107 USA  
1-801-262-3592 | [www.BioFireDefense.com](http://www.BioFireDefense.com)

*For additional information regarding our  
products and applications, please contact  
our Customer Service Department  
at BioFire Defense.*

