




**DNA Protocol  
 Basic Steps:**

- Soil
- Stool



**Preparation**

- Add sample to a SCOOP Bead Tube. See specific protocol for more information.
- Add 980 µl of **Buffer SCOOP 1A**. Add 120 µl **Buffer SCOOP 1B**.

**Lysis**

- Bead beat for 5 minutes.

**Protein Precipitation**

- Centrifuge for 2 minutes.
- Transfer liquid from bead tube to a SCOOP Receiver tube with a small transfer pipet.
- Add 250 µl **PPS** and mix by inverting 10 times.
- Centrifuge for 5 minutes.

**Bind**

- Transfer supernatant to a Mix Tube with a small transfer pipet, add 1 ml **Binding Suspension**, and invert for 2 minutes.
- Load SCOOP spin filter with mixture 4 times using Sample transfer pipet. Centrifuge each load for 2 minutes and pour off flow-through between loads.

**Wash**

- Add 500 µl **Buffer 2** to spin filter, gently stir.
- Centrifuge for 2 minutes and pour off flow-through.
- Centrifuge for 3 minutes and move spin filter to a new SCOOP Receiver tube.
- Let sit for 5 minutes with cap open.

**Elute**

- Using a Pipettor, add 200 µl of **Buffer 3** to the spin filter, gently stir.
- Centrifuge for 2 minutes.

**Keep receiver tube containing purified sample and test as soon as possible.**