

Analytical Performance of the BioFire® FilmArray® Global Fever Panel

Jared R. Helm, Corike Toxopeus, Pascal Belgique, Olivia Jackson, Alex Kelley, Cynthia D. Andjelic, Cynthia L. Phillips

BioFire Defense, LLC, Salt Lake City, UT

CONTACT INFORMATION

Jared Helm, PhD.
Jared.Helm@BioFireDefense.com
Cynthia Phillips, PhD.
Cynthia.Phillips@BioFireDefense.com

ABSTRACT

A large number of pathogens that include bacteria, viruses, and parasites can cause Acute Febrile Illness (AFI). BioFire Defense is developing the Global Fever (GF) Panel to be used on the FilmArray® System in collaboration with the U.S. Department of Defense¹ and NIAID². The FilmArray is an in vitro diagnostic test platform that combines nucleic acid purification and nested multiplex PCR for the simultaneous identification of many infectious agents in under an hour using a closed, sample-to-answer system. The FilmArray GF Panel detects and identifies nucleic acid from chikungunya virus, CCHF virus, dengue virus (serotypes 1-4), Ebolavirus, Lassa virus, Marburgvirus, West Nile virus, Yellow fever virus, Zika virus, *Bacillus anthracis*, *Francisella tularensis*, *Leptospira* spp., *Salmonella enterica* serovar Typhi and Paratyphi A, *Yersinia pestis*, *Leishmania* spp., and *Plasmodium* spp. in venous blood specimens from individuals with signs and/or symptoms of AFI or recent AFI and with known or suspected exposure to target pathogens. The LoD studies demonstrate that the GF Panel is a sensitive system that can accurately detect multiple pathogens, including Category A biothreat pathogens, and is appropriate to use in testing samples that may contain multi-analytes. Analytical specificity (exclusivity) testing demonstrates that the GF Panel is highly specific for the pathogens it is designed to detect in blood specimens. Assessment of the analytical reactivity and efficacy (inclusivity) show that the assays possess a high level of reactivity for their intended targets. A multiplex FilmArray panel could aid in rapid and actionable AFI diagnosis.


- MCS-JPEO and USAMMDA Contract No. W911QY-13-D-0080, under the NGDS program.
- NIAID Contract No. HHSN272201600002C, "Advanced Development of Multiplex Diagnostic Platforms for Infectious Diseases (Global Fever Panel)".

INTRODUCTION

The FilmArray Global Fever (GF) Panel is currently under development as a qualitative, multiplexed, nucleic acid-based test intended for use with the FilmArray 2.0 system. The FilmArray GF Panel detects and identifies bacterial, viral, and protozoan nucleic acids directly from human whole blood (EDTA) collected from individuals with signs and/or symptoms of acute febrile illness or recent acute febrile illness and with known or suspected exposure to target pathogens. The following organisms are identified using the FilmArray GF Panel: *Bacillus anthracis*, *Francisella tularensis*, *Leptospira* spp., *Salmonella enterica* serovar Paratyphi, *Salmonella enterica* serovar Typhi, *Yersinia pestis*, Chikungunya virus, Crimean-Congo hemorrhagic fever virus, Dengue virus, Ebola virus, Lassa virus, Marburg virus, West Nile virus, Yellow fever virus, Zika virus, *Leishmania* spp., and *Plasmodium* spp. (including species differentiation of *Plasmodium falciparum*, *Plasmodium vivax*, and *Plasmodium ovale*).

Figure 1. FilmArray Global Fever Pouch

The FilmArray GF pouch is a closed system disposable that stores all the necessary reagents for sample preparation, reverse transcription, polymerase chain reaction (PCR), and detection in order to isolate, amplify, and detect nucleic acid from multiple pathogens within a single clinical whole blood specimen. After sample collection, the user injects hydration solution on one side of the pouch and sample combined with sample buffer into the other side of the pouch, places the pouch into a FilmArray instrument, and starts a run. Loading the pouch takes about 2 minutes, and the entire run process takes about an hour.



LIMIT OF DETECTION

The purpose of this study is to determine the Limit of Detection (LoD) for the FilmArray GF Panel using a collection of representative organisms covering each test result. LoD₉₅ is defined as the lowest concentration of organism that can be consistently detected by the panel; analyte is detected in at least 19/20 samples (≥ 95% detected). Confirmation of the LoD is performed by testing 20 individual samples at the estimated LoD concentration.

TABLE 1. LIMIT OF DETECTION VALUES

| Organism | Strain / Source ID | Specimen | Limit of Detection (Copies/mL) | |
|-------------------------------------|---|-----------------|--------------------------------|---------|
| BACTERIA | | | | |
| <i>Bacillus anthracis</i> | Ames 35 | Live | 6.4E+01 | |
| | Ames (UCC BACI008) ² | Live | 6.4E+01 | |
| <i>Francisella tularensis</i> | Schu4 | Inactivated | 1.2E+03 | |
| | Schu4 (UCC FRAN016) ² | Live | 1.2E+01 | |
| <i>Leptospira</i> spp. | Interrogans icterohaemorrhagiae | Live | 3.9E+02 | |
| <i>Salmonella enterica</i> | Typhi | Live | 1.2E+01 | |
| | Paratyphi | Live | 6.0E+01 | |
| <i>Yersinia pestis</i> | A1122 | Live | 1.5E+02 | |
| | CO92 (UCC YERS023) ² | Live | 1.5E+02 | |
| VIRUSES | | | | |
| Chikungunya virus | R80422 | Inactivated | 5.5E+02 | |
| | B8635 (UCC Alpha031) ² | Live | 5.5E+02 | |
| | Indo23574 (UCC Alpha008) ² | Live | 5.5E+02 | |
| CCHF Virus | IbAr10200 | Inactivated | 6.4E+00 | |
| | IbAr10200 (UCC Nairo001) ² | Live | 6.4E+02 | |
| Dengue virus | DENV-1 | Hawaii | Live | 2.7E+02 |
| | DENV-2-1 | New Guinea C | Live | 3.6E+02 |
| | DENV-2-2 | Dak AR A1247 | Live | 3.6E+03 |
| | DENV-3 | H87 | Live | 1.6E+02 |
| | DENV-4 | H241 | Live | 7.6E+01 |
| Bundibugyo | 811250 (BEI) | Inactivated | 7.0E+04 | |
| | 811250 (UCC Ebola005) ² | Live | 7.0E+02 | |
| | Ivory Coast (BEI) | Inactivated | 8.3E+03 | |
| | Ivory Coast (UCC Ebola004) ² | Live | 1.8E+02 | |
| Ebolavirus | 119810 RIID (BEI) | Inactivated | 2.7E+04 | |
| | H-28(UCC Ebola003) ² | Live | 2.7E+03 | |
| | Boniface (BEI) | Inactivated | 1.1E+04 | |
| Lassa Virus | Boniface(UCC Ebola002) ² | Live | 1.1E+02 | |
| | Guéckédou/Guinea C07 (BEI) | Inactivated | 1.1E+02 | |
| | Makona(UCC Ebola027) ² | Live | 1.1E+03 | |
| Marburgvirus | Josiah (BEI) | Inactivated | 5.6E+04 | |
| | Josiah (UCC Arena002) ² | Live | 5.6E+03 | |
| | Musoke (BEI) | Inactivated | 5.0E+02 | |
| West Nile Virus | Ci67 (UCC Marbrg003) ² | Live | 5.0E+02 | |
| | Ravn (BEI) | Inactivated | 2.6E+02 | |
| | Ravn (UCC Marbrg002) ² | Live | 2.6E+02 | |
| Yellow Fever virus | NY 2001-6263 (BEI) | Inactivated | 1.6E+02 | |
| | Bz NY99 (UCC Flavi022) ² | Live | 3.2E+04 | |
| Zika virus | B-956 Uganda (BEI) | Inactivated | 1.2E+02 | |
| | 17D (BEI) | Live attenuated | 1.2E+02 | |
| <i>Leishmania donovani</i> | Asibi (UCC Flavi005) ² | Live | 1.2E+01 | |
| | PRVABC59 (BEI) | Live | 1.3E+02 | |
| PROTOZOA | | | | |
| <i>Plasmodium</i> spp. ¹ | <i>donovani</i> (BEI) | Live | 1.0E+01 | |
| | <i>falciparum</i> (clinical sample) | IPC 4884 (BEI) | Live | 1.0E+02 |
| | <i>knowlesi</i> (clinical sample) | Strain H (BEI) | gDNA | 2.2E+01 |
| | <i>malariae</i> (clinical sample) | DL517-026015 | Live | 1.9E+02 |
| <i>Plasmodium</i> ¹ | <i>vivax</i> (clinical sample) | Chesson (BEI) | Live | 7.7E+01 |
| | <i>ovale</i> (clinical sample) | N8K9QK19 | Live | 2.4E+02 |
| <i>Plasmodium</i> ¹ | <i>falciparum</i> (clinical sample) | IPC 4884 (BEI) | Live | 1.0E+02 |
| | <i>vivax</i> (clinical sample) | Chesson (BEI) | Live | 7.7E+01 |
| <i>Plasmodium</i> ¹ | <i>ovale</i> (clinical sample) | N8K9QK19 | Live | 2.4E+02 |

- Plasmodium* spp. refers to the pan assay whereas *Plasmodium* refers to the specific organism assays.
- UCC – DoD Unified Culture Collection (USAMRIID)

EXCLUSIVITY

To determine whether the FilmArray GF Panel assays cross-react with sequences from various microorganisms/viruses that may be present in clinical specimens, the analytical specificity of the panel was assessed by *in silico* analysis and by testing a broad spectrum of organisms/viruses at high concentrations. Typical stock concentrations for on-panel analytes tested are: 10⁷-10¹⁰ copies/mL for bacteria, 10⁷-10⁹ copies/mL for virus, and 10⁶-10⁸ copies/mL for protists. Both on-panel and off-panel organisms were evaluated to test inter-assay specificity and overall assay/panel specificity, respectively.

- On-panel testing consists of contrived samples spiked into sterile saline with the highest concentration of on-panel analytes that is possible based on the concentration of the organism stock (up to 10% of the total sample volume). On-panel isolates are the same as those evaluated for the LoD study (See Table 1).
- Off-panel organisms are selected based on 1) phylogenetic and/or genetic similarity to the panel analytes and assays, and 2) the possibility that the organism(s) could be present as normal flora, contaminants associated with sample collection, or pathogens in whole blood. (Organisms in red to be performed at The Center)

List of Off-Panel Organisms Tested

| BACTERIA | | |
|----------------------------------|-----------------------------------|---|
| <i>Acinetobacter baumannii</i> | <i>Enterococcus faecium</i> | <i>Rickettsia rickettsii</i> TBD |
| <i>Bacillus brevis</i> (Migula) | <i>Francisella</i> (4 strains) | <i>Rickettsia typhi</i> TBD |
| <i>Bacillus</i> (12 strains) | <i>Haemophilus</i> | <i>Salmonella bongori</i> |
| <i>Bacteroides fragilis</i> | <i>Klebsiella oxytoca</i> | <i>Salmonella enterica</i> (18 Strains) |
| <i>Bordetella bronchiseptica</i> | <i>Legionella pneumophila</i> | serovar Dublin (Paratyphi) |
| <i>Borrelia burgdorferi</i> | <i>Leptospira</i> (6 strains) | serovar Manchester (Typhi) |
| <i>Brucella melitensis</i> | <i>Listeria monocytogenes</i> | <i>Serratia marcescens</i> |
| <i>Burkholderia</i> (3 strains) | <i>Mycobacterium tuberculosis</i> | <i>Staphylococcus aureus</i> |
| <i>Chlamydia pneumoniae</i> | <i>Mycoplasma pneumoniae</i> | <i>Streptococcus</i> (3 strains) |
| <i>Chlamydia psittaci</i> TBD | <i>Neisseria meningitidis</i> | <i>Treponema pallidum pallidum</i> |
| <i>Clostridium</i> (4 strains) | <i>Proteus mirabilis</i> | <i>Vibrio cholerae</i> |
| <i>Coxiella burnetii</i> | <i>Pseudomonas aeruginosa</i> | <i>Yersinia</i> (11 Strains) |
| <i>Enterobacter aerogenes</i> | <i>Rickettsia prowazekii</i> TBD | |
| <i>Enterococcus faecalis</i> | | |

| VIRUSES | | |
|-----------------------------------|---|---|
| Adenovirus 1, 3, 5 | Human Herpesvirus 6B | Powassan virus TBD |
| Aura virus | Human Immunodeficiency Virus Type 1, 2 | Rabies virus |
| Avalon virus (in silico) | Human T-Lymphotropic Virus Type 1, 2 | Rift Valley Fever Virus |
| Barnah Forest virus | Influenza A H1N1-2009 | Ross River virus |
| Bas-Congo (in silico) | Influenza A H3N2 | Human respiratory syncytial virus |
| Bunyamwera virus | Influenza B virus | Rubella virus |
| Coronavirus | Japanese encephalitis virus | Saint Louis encephalitis virus |
| Cytomegalovirus | Junin virus | Sabia virus (in silico) |
| Dugbe virus TBD | Lymphocytic choriomeningitis virus | Semliki Forest virus TBD |
| Eastern equine encephalitis virus | Mumps virus | Sindbis virus |
| Enterovirus, HEV-71 | Measles virus | Spondweni virus |
| Epstein Barr virus | Metapneumovirus | Tickborne encephalitis virus |
| Flexal virus | Middelburg virus TBD | Una Virus |
| Guanarito virus | Mopeia virus TBD | Usutu virus |
| Hantaan virus | Murray Valley encephalitis virus TBD | Vaccinia Virus |
| Hazara virus | Ormsk hemorrhagic fever | Variella Zoster virus |
| Hendra Virus | O'Nyong Nyong virus (CHIKV) | Variola majora (In Silico) |
| Hepatitis A, B, C virus | Parvovirus | Venezuelan Equine Encephalomyelitis virus |
| Herpes Simplex Virus Type 2 | Pirital virus (In Silico) | Western Equine Encephalomyelitis virus |
| HPIV-1, 3 | | |
| Hughes virus | | |

| FUNGI AND PROTISTS | | |
|--|--|---|
| <i>Babesia microti</i> | <i>P. cynomolgi</i> (plus vivax/ovale) | <i>Toxoplasma gondii</i> |
| <i>Candida fasciculata</i> (Leish) | <i>P. fieldi</i> (plus vivax/ovale) | <i>Trypanosoma</i> (3 strains) |
| <i>Cyclospora cayentanensis</i> | <i>P. fragile</i> (plus vivax/ovale) | <i>Aspergillus fumigatus</i> |
| <i>Leptomonas seymouri</i> (Leish) | <i>P. inui</i> (plus vivax/ovale) | <i>Cryptococcus neoformans</i> var. <i>grubii</i> |
| <i>P. berghii</i> (plus vivax/ovale) | <i>P. simiovale</i> (plus vivax/ovale) | |
| <i>P. brasilianum</i> (plus vivax/ovale) | <i>Schistosoma mansoni</i> | |

INCLUSIVITY

To ensure the FilmArray GF Panel is inclusive for the genetic variation expected for each analyte, the analytical reactivity of the assays was evaluated by testing multiple isolates per analyte. Isolates are selected based on the availability of live and inactivated stocks, genetic, temporal and geographic diversity, and clinical relevance of the various species, strains, subspecies, serotypes, genotypes and genetic variants available for testing. In addition, *in silico* data (sequence searches and alignments to assay primers) are also used to support the inclusivity of the FilmArray GF Panel assays.

Samples are prepared by spiking analytes into whole blood from healthy donors (obtained from Bioreclamation IVT repository) at a concentration near 3x LoD.

TABLE 2. INCLUSIVITY

| FilmArray Global Fever Panel Analyte | # Isolates Detected/ Isolates Tested |
|--|--------------------------------------|
| BACTERIA | |
| <i>Bacillus anthracis</i> | 4/4 |
| <i>Francisella tularensis</i> | 5/5 |
| <i>Leptospira</i> spp. | 16/16 |
| <i>Salmonella enterica</i> serovar Paratyphi A | 3/4 |
| <i>Salmonella enterica</i> serovar Typhi | 5/5 |
| VIRUSES | |
| Chikungunya virus | 3/3 |
| DENV1 | 7/7 |
| DENV2 | 6/7 |
| DENV3 | 5/5 |
| DENV4 | 7/7 |
| Ebolavirus | Zaire 2/2 |
| Marburgvirus | 3/3 |
| West Nile virus | 3/3 |
| Zika virus | 2/2 |
| Zika virus | 7/7 |
| PROTOZOA | |
| <i>Leishmania</i> spp. | 9/9 |
| <i>Plasmodium</i> spp. | 17/17 |

REPRODUCIBILITY

This study evaluated the reproducibility of the FilmArray GF Panel test results. The study was designed to evaluate the potential for run-to-run or day-to-day variation in panel performance, as well as potential variation linked to test location (site), reagent lot, operator, or instrument for samples containing analytes at and near their LoD.

Three contrived whole blood samples were spiked with different mixtures of six representative panel analytes, two bacteria (*Leptospira interrogans* and *Salmonella enterica enterica* serovar Typhi), two viruses (dengue virus and Zika virus), and two protozoa (*Leishmania donovani* and *Plasmodium falciparum*). One sample was spiked at a Moderate Positive (3x LoD) level, and another sample at a Low Positive (1x LoD) level; the third sample was not spiked. The three samples were aliquoted and distributed between three test locations. Six replicates of each sample were tested at each location on five different days, with two operators and two instruments, providing a total of 90 replicate test results per sample. In total, valid results were obtained from testing 270 sample aliquots to demonstrate the reproducibility of the FilmArray GF Panel.

The overall test results for six analytes showed a 97.5% agreement of the observed test results with the expected Detected or Not Detected results (1843/1890), with a 95% confidence interval of 96.7-98.1% (Table 3).

Table 3. Results for Reproducibility testing on the Global Fever Panel

| Moderate Positive (3X) | 626/630 (99.4%) | Overall Agreement (All analytes, all levels) | 1843/1890 (97.5%) [96.7-98.1%] |
|------------------------|-----------------|--|--------------------------------|
| Low Positive (1X) | 587/630 (93.2%) | | |
| Negative | 630/630 (100%) | | |

INTERFERING SUBSTANCES

This study was designed to evaluate the potential for certain substances to interfere with the performance and accuracy of FilmArray GF Panel test results. These substances included; (i) endogenous substances, (ii) exogenous substances, (iii) microorganisms such as viral and bacterial species, and (iv) technique specific substances that may be introduced into a sample (intentionally or not) during routine laboratory handling. Interference was not detected under any conditions tested (Table 4).

Table 4. Results for Potentially Interfering Substances Tested on the Global Fever Panel.

| Potentially Interfering Substances | | | | | | |
|------------------------------------|----------------------------|---------------------------------|-----------------------------|------------------------------------|-------------------------------------|-----------------|
| EXOGENOUS SUBSTANCES (EX) | ENDOGENOUS SUBSTANCES (EN) | COMPETITIVE MICROORGANISMS (MO) | BLOOD COLLECTION TUBES (BT) | TECHNIQUE SPECIFIC SUBSTANCES (TS) | | |
| Artemether Lumefantrine | Clindamycin | Osetamivir | Albumin | Corynebacterium diphtheriae | Citrate (sodium) | Bleach |
| Acetaminophen | Cortisol | Pentamidine | Bilirubin (Conjugated) | Cytomegalovirus (CMV) AD-169 | Acid-citrate-dextrose (ACD) | Acetone |
| Amoxicillin | Cycloserine | Prednisolone | Bilirubin (Unconjugated) | Epstein-Barr virus | EDTA in excess (5x) | Chloroform |
| Amphotericin B | Doxycycline | Prednisone | Cholesterol (total) | Escherichia coli | Heparin | DMSO |
| Artesunate | Fluconazole | Proguanil | Glucose | Haemophilus influenzae | Sodium polyanethane-sulfonate (SPS) | Ethanol |
| Aspirin (Acetylsalicylic Acid) | Gentamicin | Ribavirin | Hemoglobin | Herpes Simplex virus | | Methanol |
| Atovaquone | Ibuprofen | Sulfamethoxazole | Immunoglobulins | Human Immunodeficiency virus (HIV) | | Povidone-iodine |
| Azithromycin | Isoniazid | Tenofovir | Triglycerides | Klebsiella pneumoniae | | Saline |
| Ceftriaxone | Mefloquine | Vancomycin | White Blood Cells | Plasmodium vivax | | TRizol |
| Ciprofloxacin | Meropenem | | | Staphylococcus epidermidis | | |

SAMPLE STORAGE AND TRANSPORT

The purpose of this study was to validate that accurate FilmArray GF Panel test results can be obtained when testing human whole blood (WB) specimens that have been stored at ambient room temperature for one day, at 4°C for up to seven days, or in an ultra-low temperature freezer. The aim of the study was not to establish the limits of analyte stability and sample storage, but to confirm that analytes are stable and detected by the panel when stored under the selected conditions.

Testing was conducted using samples composed of human WB, each contrived with three representative GF Panel analytes at a concentration equal to 3x their LoD. To evaluate the impact of storage, 10 replicates were tested at each storage condition.

All analytes were Detected in ≥9/10 replicates at each storage condition. The robust detection of representative mixes observed in this study supports the conclusion that accurate FilmArray GF Panel test results can be obtained from human WB specimens stored in the conditions outlined.

SUMMARY

The development of a multiplex FilmArray panel would aid in rapid and actionable AFI diagnosis. The FilmArray GF Panel provides a broad spectrum analysis of target pathogens in a sensitive and specific manner.

- LoD values show sensitivity levels at or below clinically relevant concentrations.
- On- and off-panel exclusivity show very limited evidence of cross-reactivity.
- Inclusivity maintains expected coverage of relevant target pathogens.
- Reproducibility study ensures consistent test results.
- No interference from substances tested show panel robustness.
- The GF Panel provides accurate results for Whole Blood samples stored under several standard conditions.

Data presented are from assays that are Investigational Use Only (IUO) and have not been cleared or approved for diagnostic use.